

Globules

Quality of Globules:

Deeplaid Laboratories Ltd. using the pure sugar of milk for its manufacturing i.e coating of sugar cubes, whereas other brands might be using normal sugar powder. Why we say this is that Deeplaid Laboratories globules do not dissolve. The integrity is well maintained. When you add the dilution to the globules, the dilution travels immediately to the bottom of the bottle, making the globules moistened homogeneously all through with very few drops as 1-3 drops for a 1 drum vial containing Deeplaid Laboratories globules. In case of other brands, as you add the dilution, you will find the globules start dissolving and the dilution do not reach to the bottom of the vial moistening the globules, and if you add more drops of dilution, you will see the globules dissolve forming a paste. So, that's how Deeplaid Laboratories globules follow strictly the Hahnemanian principles of minimum dosage.

Quality Maintain Procedure:

- Not at all touched by hand at any stage of production.
- Homogeneity variation of size is very minimal.
- Better absorption.
- Test solvent reaches even the centre of the globules during quality control and quality assurance tests.
- Entire process is done under highly hygienic condition.
- Only Pharma grade sugar is used .
- Raw materials are thoroughly checked and ensured for highest quality.
- Controlled Temperature Treatment Technique (CTTT) ensures highest softness of globules Colour of the globule is pure natural white.
- No external material is added to artificially make the globules brighter.

Size: Pill -20, 30, 40

Color: White

Presentation: 375gm in poly bag.